

# OPHTHALMIC GEL



eyeDRO IS A TRIPHASIC POLYMERIC GEL DESIGNED TO PROTECT AND HYDRATE THE CORNEAL SURFACE DURING CATARACT AND VITREORETINAL SURGERY

eyeDRO REPLACES MULTIPLE BSS OR VISCOELASTIC APPLICATIONS IN CATARACT AND VITREORETINAL SURGERY

## TRIPHASIC POLYMERIC GEL

REF.	PRODUCT	VOL.	QTY
EDO 001	eyeDRO advanced polysaccharide blend	2 ml	10



EDO 001

# OPHTHALMIC GASES



## MULTIDOSE CYLINDER

↑ MORE GAS

↓ reduced dimensions

REF.	PRODUCT	VOL.	QTY
GOT 007-00	GOT Multi SF <sub>6</sub> pure sulfur hexafluoride	75 ml	1
GOT 008-00	GOT Multi C <sub>2</sub> F <sub>6</sub> pure hexafluoroethane	75 ml	1
GOT 009-00	GOT Multi C <sub>3</sub> F <sub>8</sub> pure octafluoropropane	75 ml	1

Each pack contains 1 gas cylinder + 1 gas injection set  
+ 1 needle and 1 patient's wristband

## SET FOR REPEAT USE OF OPHTHALMIC GAS

REF.	PRODUCT	QTY
GIS 003-00	GIS gas injection set made of filter, needle, and patient's wristband	25



GOT 007-00



GOT 008-00



GOT 009-00



GIS 003-00

# PERFLUOROCARBONS

## SYRINGE

REF.	PRODUCT	VOL.	QTY
HPF 019-00	HPF8 high purity perfluoro-n-octane	5 ml	1
HPF 020-00	HPF8 high purity perfluoro-n-octane	7 ml	1
HPF 021-00	HPF10 high purity perfluorodecalin	5 ml	1
HPF 022-00	HPF10 high purity perfluorodecalin	7 ml	1



HPF 019-00



HPF 020-00




HPF 001-00



HPF 003-00

## VIAL

REF.	PRODUCT	VOL.	QTY
HPF 001-00	HPF8 high purity perfluoro-n-octane	5 ml	1
HPF 002-00	HPF8 high purity perfluoro-n-octane	7 ml	1
HPF 003-00	HPF10 high purity perfluorodecalin	5 ml	1
HPF 004-00	HPF10 high purity perfluorodecalin	7 ml	1



# SILICONE OILS



## SYRINGES

REF.	PRODUCT	VOL.	QTY
RSO 006-00	RS-OIL ECS silicone oil 1000 cS	10 ml	1
RSO 007-00	RS-OIL ECS silicone oil 5000 cS	10 ml	1
RSO 010-00	RS-OIL ECS silicone oil 1300 cS	10 ml	1
RSO 011-00	RS-OIL ECS silicone oil 5700 cS	10 ml	1
RSO 013-00	RS-OIL ECS silicone oil 2000 cS	10 ml	1

## VIAL

REF.	PRODUCT	VOL.	QTY
RSO 002-00	RS-OIL silicone oil 1000 cS	15 ml	1
RSO 003-00	RS-OIL silicone oil 1000 cS	10 ml	1
RSO 004-00	RS-OIL silicone oil 5000 cS	15 ml	1
RSO 005-00	RS-OIL silicone oil 5000 cS	10 ml	1
RSO 008-00	RS-OIL silicone oil 1300 cS	15 ml	1
RSO 009-00	RS-OIL silicone oil 5700 cS	15 ml	1
RSO 012-00	RS-OIL silicone oil 2000 cS	15 ml	1



RSO 006-00



RSO 007-00



RSO 010-00



RSO 011-00



RSO 013-00



RSO 002-00



RSO 004-00



RSO 008-00



RSO 009-00



RSO 012-00